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DEPARTMENT OF INTERNAL MEDICINE
SEVERANCE HOSPITAL
YONSEI UNIVERSITY COLLEGE OF MEDICINE
SEOUL, 120-752, KOREA

250 Seongsanejo, Seodaemun-Ku,
(134 Shinchon-Dong) Seoul 120-752 Korea
TEL : (02)-2-2228-1930, 1931
FAX : (02)-2-393-6884
E-mail : medi@yuhs.ac

May, 2nd, 2011

1. Test: Blocking effects of allergen penetration of NOSK®
2. Summary:

We certify that the intranasal filter, NOSK® can effectively block the penetration of house dust mite (84.3% reduction of house dust mite protein by protein assay method, 73% reduction of house dust mite allergen by ELISA inhibition assay) and oak pollen (95.0% reduction of oak pollen allergen by protein assay method) particles as demonstrated through using the customized airflow chamber.

Prof. Jung-Won Park, MD., PhD.

Dept. Internal Medicine, Yonsei University, College of Medicine

Seoul, Korea 120-752

Tel; 82-2-2228-1961

Fax; 82-2-2227-7932

Email; parkjw@yuhs.ac

Test Report

Test object: Intranasal filter, NOSK®

Test: Blocking effects of allergen penetration via intranasal filter, NOSK®

Materials and methods

I. Materials

Two kinds of allergenic particle were used. One is house dust mite culture media powder (HDMCM, culture media of *Dermatophagoides farinae*) that containing abundant house dust mite allergens was kindly provided by the Arthropodes of Medical Importance Resource Bank & Institute of Tropical Medicine. The second is oak tree pollen powder provided by Institute of Allergy in Yonsei University College of Medicine.

We made a customized chamber for this test. This chamber consisted of two compartments which were divided by a 72-holed acryl plate (Fig 1). We placed 36 pairs of NOSK® on the 72-holes at acryl plate. An anemometer was placed in front compartment to measure the speed of airflow. A filter paper was placed at opposite side of the rear chamber to measure the amount of penetrated allergens (Fig. 2).

II. Methods

Test for blocking effects of allergen penetration via NOSK®

1) Test procedure

For each allergenic powder, tests were conducted two times with/without NOSK® at 72-holed at acryl plate. Initially 20 g of allergenic powder was loaded at the front chamber without NOSK®. Once turned on the air blower, continuous airflow was maintained with 0.53~1.92 m/sec speed for 30 minutes. The speed of airflow of the chamber was changed as the nasal filters were fitted (0.84~1.92 m/sec without NOSK®, 0.53~0.72 m/sec with NOSK®). Filter paper at rear chamber was used to collect the penetrated particles. We extracted proteins from the filter paper and used following assays.

2) Extraction of allergens and protein assay

Each filter paper from rear chamber was incubated with 2 mL of 0.1M bicarbonate buffer solution (pH 8.0) in 50 mL tube using orbital shaker. Protein extraction was performed for 18

hours at room temperature. Finally we collected supernatants after centrifugation of the extracts (3000 rpm, 20 min) and use them in following assays. We performed Bradford protein assay (Bio-rad, CA) to measure the total protein content in the extract solution.

3) ELISA inhibition assay

96-well ELISA plate was used. We coated each well with 100 μ L of *D. farinae* crude extract in carbonate buffer (pH 9.6) (10 μ g/mL concentration) and incubated overnight at 4°C. Pooled serum of *D. farinae* sensitized patients (n=5) were incubated with serially diluted protein extract solution from the filter papers as inhibitor for 24 hours. Blocking with bovine serum albumin, we loaded inhibitor mixed sera 100 μ L in each well as primary antibody. One hour later, 1:1000 diluted anti-human IgE antibody (Vector, CA) was added and visualized with tetramethylbenzidine and microwell peroxidase substrate (KPL, USA). Optical density was read by spectrophotometer (Molecular devices, Korea) with 450 nm UV light. Finally we calculate the % of inhibition according to the added inhibitor concentration.

III. Results

1) Gross inspection

After the procedure we observed the inside of both chambers. In the front chamber we could easily find allergenic particles those could not pass through NOSK® filter (Fig. 3). On the other chamber, it was hard to find allergenic powder in gross inspection.

2) Bradford protein assay for the extractions

The protein content of house dust mites culture media in the extraction was decreased from 1,200 mg/ml to 189 mg/ml (84.3% reduction) by adapting NOSK® filters. The oak pollen protein content was decreased from 1,200 mg/ml to 60 mg/ml (95.0% reduction) (Fig. 4).

3) ELISA inhibition assay for the extractions

We observed dose-dependent increment of inhibition ratio in both settings with/without NOSK® filter. Especially at 1:1 concentration between atopic pooled serum and the extract solution from the filters with NOSK® showed 38% inhibition of specific IgE binding to *D. farinae* allergens. Its potency was same as 1:0.27 concentration of none filter condition. This result suggested that NOSK® filter could reduce penetration of house dust mite allergen by 73% compared by non-filter condition.

IV. Figures


Figure 1. 72-holed acrylic plate that divide chamber to two compartments (36 pairs of NOSK® filters were placed.)


Figure 2. Customized chamber used in this experiment. An anemometer and a filter paper were placed.


Figure 3. Gross inspection of 72-holed acrylic plate. Large amount of allergenic particles could not pass through NOSK® filters were noted.


Figure 4. Total protein amounts measured by Bradford assay were remarkably decreased from the baseline in the experiments with HDMCM and oak pollen (84.3% and 95.0%), respectively. *HDMCM; house dust mite culture media particle

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